


18, rue Marie et Pierre Curie
92260 FONTENAY AUX ROSES
Phone : 01 46 60 61 51
fax : 01 78 76 74 18
http : // www.fime.fr



MANUAL READER FOR USE IN T.L. DOSIMETRY



- Windows executable software
- Seven decade range
- Save all files and measurement results
- Automatic assessment of statistical values (mean, standard deviation)
- Automatic calibration of the acquisition line using a high accuracy internal reference light ($< 0.2\%$)
- Glow curve printing
- Interchangeable heating planchette for all types of dosimeters (powder pellets, microrods, Teflon disc, etc..)
- Interchangeable optical filters
- Maximum temperature : $500\text{ }^{\circ}\text{C}$